

Pilot investigation into effect of an antimicrobial dressing on the cell migration effects of decellularized porcine placental extracellular matrix in viable wounded human skin *ex vivo*

Alicia Hughes BSc¹, Leyla Zilic PhD¹, Rocio Burgos-Amador PharmD¹, Daniel Metcalf PhD¹

¹Advanced Wound Care R&D, Convatec Ltd, Deeside, UK

Introduction

- Decellularized extracellular matrix (ECM) medical devices can facilitate healing of hard-to-heal wounds¹
- Porcine placental ECM is an ECM scaffold containing collagen, fibronectin, laminin, elastin, hyaluronic acid, and glycosaminoglycans, while being largely free of cells, cell debris, and DNA²
- However, as hard-to-heal wounds are often compromised by microbial bioburden, so the use of an antimicrobial dressing over ECM products may be desired
- An advanced antimicrobial dressing* designed to manage wound bioburden, including surface-associated and aggregated microorganisms, has been shown to facilitate healing in hard-to-heal wounds³, and as the dressing component in Wound Hygiene⁴

STUDY OBJECTIVE

To examine the effects of a novel porcine placental decellularized ECM product* covered with an advanced antimicrobial wound dressing† on cell migration in a wounded human skin in a viable *ex vivo* model in a pilot study using immunolabeling techniques

Methods

Test devices

- Porcine placental ECM*
- Non-adherent wound contact layer†
- Carboxymethylcellulose dressing with ionic silver, surfactant, and chelator (antimicrobial dressing)§

Protocol

- Human *ex vivo* wounded skin models were treated with the ECM product alone, the antimicrobial dressing alone (fully hydrated with simulated wound fluid), or a combined construct (Figure 1) containing the ECM product, a non-adhesive wound contact layer†, and the hydrated antimicrobial dressing)§
- Untreated control, antimicrobial dressing-treated, porcine ECM-treated, and construct-treated wounded skin samples were incubated for up to 6 days at 37°C. Skin samples were fixed and embedded before cross sectioning
- Histological slides were stained with Cytokeratin 17 primary antibody (1:200). Cytokeratin 17 was selected as a marker for wound healing as it is associated with cellular migration
- Primary antibody was removed, and fluorescently labelled with the secondary antibody, Alexa Fluor 647 (1:500) ($\lambda_{ex} = 650$, $\lambda_{em} = 671$)
- Additionally, the slide was counterstained with DAPI (1:1000) ($\lambda_{ex} = 350$, $\lambda_{em} = 465$) to visualize cell nuclei. Slides were observed under a LSM800 Zeiss Confocal Laser Scanning Microscope

Results

- Compared to the untreated control (Figure 2), Cytokeratin 17 was still observed in the antimicrobial dressing-treated wounded skin sample (Figure 3)
- Notably, increased upregulation of Cytokeratin 17 was observed in the porcine ECM-treated (Figure 4) and the construct-treated wounded skin samples (Figure 5)




Figure 2. Untreated control wounded skin




Figure 3. Antimicrobial dressing-treated skin




Figure 4. Porcine ECM-treated wounded skin




Figure 5. Construct-treated wounded skin

Discussion

- This *ex vivo* pilot study suggests that the wound healing potential of the porcine placental decellularized ECM device is not impacted by the presence of an advanced antimicrobial dressing
- Furthermore, the antimicrobial and physical properties of a silver-containing antimicrobial dressing, when applied over a non-adherent layer, does not affect the ability of a porcine placental decellularized ECM device to encourage cell migration in a wounded human skin model
- Further studies in similarly complex models using a range of microscopic and immunological techniques, such as Ki67 (an antibody associated with cell proliferation), may help to confirm these pilot findings

References & Footnotes

1. Cramer M, Badyal, SF. Extracellular Matrix-Based Biomaterials and their influence upon cell behavior. *Ann Biomed Eng* 2020; 48: 2132-2153.
2. FDA 510(k) summary K193552: www.accessdata.fda.gov/cdrh_docs/pdf21/K211902.pdf (accessed March 2024).
3. More Than Silver™ Made Easy. *Wounds Intl* 2019 (available to download at www.woundsinternational.com [accessed March 2024]).
4. Murphy et al. International consensus document. Defying hard-to-heal wounds with an early antibiofilm intervention strategy: wound hygiene. *J Wound Care* 2020; 29(Suppl 3b): S1-28.

CONCLUSION

These findings suggest that the wound healing potential of the porcine placental decellularized ECM device is not impacted by the presence of an advanced antimicrobial dressing




Figure 1. ECM-non-adherent antimicrobial dressing construct-treated human wounded skin model

*InnovaMatrix® AC (Convatec Inc)
†Mepitel (Molnlycke)
§Aquacel® Ag Advantage dressing (Convatec Inc)
Editorial support was provided by Kenny Tran (Convatec Ltd.)