

Profiling of cyclin-dependent kinase (CDK) inhibitors: use case of a novel DIANA method for efficient kinase inhibitor discovery



Ilona Kalasova, Lucia Bialesova, Alena Uvizl, Pavel Svec, Jiri Schimer, Vaclav Navratil, Martin Dienstbier, Jan Tykvart

+/-

+

Overview

- **DIANA** (DNA-linked A novel Inhibitor ANtibody Assaytechnology enabling efficient inhibitor profiling and HTS (1).
- Convenient determination of enzyme inhibition constants (K_i) in a 384-well plate format.
- Reliable determination of K_i from a single well.



Introduction

Cyclin Dependent Kinases (CDKs) are serine/threonine kinases, members oft he CMGC kinase family, whose activity is regulated by their binding partners, cyclins (2). To this day, 20 CDKs have been identified (3). They play a role in essential cellular processes such as regulation of the cell cycle, transcription, and pre-mRNA splicing.

Here, we present DIANA CDK panel – a set of microplatebased binding assays for CDKs based on the patented DIANA (DNA-linked inhibitor antibody assay) technology. These assays can be used for convenient measurement of inhibitors' K_i and selectivity profiling of target compounds against the whole family of CDKs. Alternatively, DIANA CDK assays can be applied as robust HTS (high-throughput screening) method а to search for new active compounds against CDKs.

Selectivity profiling of target compounds against the family of CDKs.

1 | Solid phase binding 2 | Probe + Inhibitor competition/binding 3 | qPCR readout \sim \checkmark K_i calculated **Cost effective Highly sensitive** from a single well Test compound 1. Antibody activity **DIANA probe** Reporter DNA Linker 🔨 2. Enzyme Highly active compound Inhibitor Cq Active compound Detected Test probe compound activity binding measurement

Features of DIANA CDK panel

- Based on the proprietary patented DIANA technology.
- 384-well plate format → **high-throughput.**
- Low consumption of reagents → **cost-effective.**
- Quantitative PCR-based readout → K_i from single well.
- Any CDK included in the DIANA CDK inhibitor

profiling panel can be used in **HTS campaigns** (Z' > 0.7) using our in-house compound libraries.

Parameters of DIANA CDK panel

CDK	Dynamic range (fold change)	MIN K _i (nM; 1 µM cpd used)	MAX K _i (nM; 1 µM cpd used)
CDK1/cycB1	15	110	1700
CDK2/cycA2	2100	1.0	2200
CDK3/cycE1	1600	1.4	2300
CDK4/cycD1	7.7	210	1600
CDK4/cycD2	13	130	1600
CDK4/cycD3	45	44	1900
CDK5/p25	550	3.4	1900
CDK5/p35	700	2.6	1800
CDK6/cycD1	17	77	1300
CDK6/cycD2	5.4	330	1700
CDK6/cycD3	26	77	2000
CDK7/ cyc H1	15	100	1600
CDK7/ cyc H1/ MNAT1	1500	0.90	1300
CDK8/cycC	13	93	1200
CDK9/cycT1	320	2.8	910
CDK9/cycT2	290	5.3	1500
CDK10/cycQ	75	26.0	1900
CDK11A	49	25	1200
CDK12/cycK	62	35	2200
CDK13/cycK	61	27	1700
CDK14/cycY	330	5.7	1900
CDK16/cycY	35	58	2000
CDK17/cycY	57	39	2200
CDK18/cycY	7.5	260	1900
CDK19/cycC	160	12	1900

Methods

Cell cycle CDKs 📃 Transcriptional CDKs

Contact information

- 0 Průmyslová 596, 252 50 Vestec, Czech Republic
- Tel: +420 212 247 340

e-mail: info@dianabiotech.com

(in)

DIANA Biotechnologies

Contact us at: www.dianabiotech.com

panel in your drug

discovery pipeline?

Conclusions ·

DIANA method enables robust, precise, and cost-effective inhibitor profiling.

- DIANA CDK panel comprises 25 CDK/ cyclin complexes.
 - Results of the DIANA CDK panel were validated with two reference methods.
 - K_i determination from a single well.

References

- 1. Navrátil et al., Nucleic acids research 2017
- 2. Malumbres et al., Genome Biology 2014
- 3. Lukasik et al., Int. J.Mol. Sci. 2021