From Microcavities to Morphological Insights: High-Throughput 3D Organoid Imaging and Analysis with Gri3D® and CQ1 Ballabio, M.¹, Wieland, F. A.², Clapés Cabrer, M.¹, Meyer, M.³, Hoehnel-Ka, S.¹ ¹SUN bioscience SA, Lausanne, Switzerland - ²Yokogawa Life Science Europe, Basel, Switzerland - ³Doppl SA, Lausanne, Switzerland

INTRODUCTION

Organoids bridge traditional cell culture and animal models, offering promise in various fields, and are traditionally cultured in solidified extracellular matrix (ECM). This method presents multiple challenges, including high-content imaging in 3D due to sample heterogeneity and variations in focal planes [1].

ECM-dome 3D culture

* U-bottom microcavities



Feature	ECM-dome	Gri3D®
HOMOGENEITY		
LOCATION (x,y,z)		
AUTOMATION		
CELL-CELL INTERACTIONS		

To tackle these hurdles, we introduce Gri3D[®]: an innovative hydrogel microwell 96WP system designed for uniform cell seeding, efficient aggregation, and the generation of individual microtissues in suspension-like conditions [2]. Importantly, the resulting microtissues are tesselated all within the same focal plane, enabling simultaneous high-resolution imaging. Here, we combined the plate with Yokogawa's CQ1, a dual-microlens spinning disk imaging confocal microscope, which has proven to be successful in high-content imaging of 3D spheroid models [3]. In this study, we demonstrate high-content imaging of mouse intestinal organoids cultured in Gri3D[®] plates using Yokogawa's CQ1 high-content imaging system and CellPathfinder software, facilitating precise segmentation and analysis.



METHODS

Mouse intestinal organoids are cultured in Gri3D[®] plates for five days. Subsequently, the microtissues are fixed with paraformaldehyde, stained, and imaged directly on an imagingbottom plate, with Yokogawa's CQ1. The transparent hydrogel enables high-resolution imaging of organoids in microcavities, with Z-stack imaging revealing intricate 3D structures. CellPathfinder software facilitates segmentation and analysis, yielding insights into organoid characteristics such as nuclei count, epithelial thickness, and lumen definition.

Model: Mouse Intestinal organoids

<u>Gri3D® Format</u>: 400 µm, 121 µwells - imaging bottom - black walls Seeding density: ~200 cells per microwell (24.2k cells per well) Cell types: Mouse intestinal stem cells (mISC) <u>Culture time:</u> 5 days





RESULTS



Figure 1) Organoid fluorescence imaging and segmentation of one single Gri3D® well, containing 121 3D microtissues. Left: 4x stitched images taken on a Yokogawa CQ1. Maximum intensity projection of stained mouse intestinal organoids. DAPI in Blue and Phalloidin in Magenta. Right: Segmentation, labeling, and counting of the organoids in Yokogawa's CellPathfinder software. Scale bars 400 µm.

Gri3D[®] 96WP 400 µm microwells are used to scale up the generation of intestinal organoids. 121 microtissues are formed within one well, all in the same focal plane. After five days, the organoids are incubated with the antibodies of interest, fixed, and immunostained on Gri3D® following an automatable protocol (Figure 1. Left). The images are analyzed with CellPathfinder, allowing to automatically count,



label, and segment the 121 microtissues (Figure 1. Right).

Because the hydrogel is transparent (95% water), high-magnification (20x) Z-stack imaging can be performed directly on-plate using the CQ1, and further 3D analyses can be conducted with CellPathfinder to extract numerous features. In this instance, we successfully segment the cell nuclei, measure the epithelium thickness, and accurately define the hollow lumens inside the mouse intestinal organoids in 3D (Figure 2).



The integration of Yokogawa's CQ1 high-content system with SUN bioscience's Gri3D[®] plates represents a significant advancement in **3D imaging for organoid research**. The combined use of Yokogawa's cutting-edge imaging technology and SUN bioscience's innovative Gri3D[®] plates addresses previously encountered challenges by enabling faster image acquisition and more accurate data. This synergistic approach opens new avenues for advancing organoid-based research, paving the way for deeper insights into complex biological systems and accelerating the development of novel therapeutics.

REFERENCES	ACKNOWLEDGEMENTS
 Fatehullah, A., Tan, S. H. & Barker, N., 2016. Organoids as an in vitro model of human development and disease. <i>Nature Cell Biology</i> 18, 246–254. Brandenberg, N., Hoehnel, S., Kuttler, F., et al., 2020, 'High-throughput automated organoid culture via stem-cell aggregation in microcavity arrays', <i>Nat Biomed Eng</i> 4, 863–874. Wardwell-Swanson J, Suzuki M, Dowell KG, et al., 2020. 'A Framework for Optimizing High-Content Imaging of 3D Models for Drug Discovery'. <i>SLAS Discovery</i> 25(7):709-722 	The organoids were grown, fixed, and stained in the Doppl SA laboratories by Ballabio M., and Clapés M. The images were acquired at the Yokogawa facilities in Basel, by Wieland F.A., Ballabio M., Clapés M., and Meyer M. Gri3D plates were provided by SUN bioscience. Ballabio M., Clapés M., and Hoehnel-Ka S. conceptualized and created this poster and all rendered visuals.
SUNBICSCIENCE YC	NOGANA DOCL