# A common laboratory challenge: Harmonizing the HLA-DPB1 antibody specificity nomenclature, and UNOS requirements, for unacceptable antigens in HistoTrac<sup>TM</sup> using software automation

Phillip R. Fisher<sup>1, 2</sup>, Neng Yu<sup>3</sup>, Kesherwani, Varun<sup>3</sup>, Efi Melista<sup>1</sup>, Patrick Ching<sup>1</sup>, Susie Brown<sup>1</sup> <sup>1</sup>Transplant Diagnostics Division, Thermo Fisher Scientific, Statesville, NC, United States. <sup>2</sup>ASHI, Statesville, NC, United States. <sup>3</sup>HLA Services, American National Red Cross, Dedham, MA, United States.

### Introduction

The HLA laboratories in the US supporting solid organ transplantation use an established process to upload unacceptable antigen (UA) results to UNET following ASHI & OPTN(UNOS) Standards. UNOS requires the DPB1 antibody specificities to be reported at a two-field allele level. This is in part to meet current typing methodologies employed by HLA labs as well as to ensure patient safety, organ allocation, and accurate reporting of DPB1 antibodies that can affect graft survival. The current HLA antibody single antigen bead (SAB) assays offered by OneLambda<sup>™</sup> uses HLA Fusion<sup>™</sup> software for analysis, which uses a terminology different from what is acceptable under current standards, forcing the lab staff to manually edit every DPB1 antibody assignment. Since the terminology used is outdated for todays standards, it could result in inappropriate organ allocation, unexpected positive crossmatches, increased cold ischemia time, and the potential for transplanting incompatible organs. This work aims to develop a solution that harmonizes the DPB1 antibody assignment nomenclature in HistoTrac to meet ASHI and UNOS requirements automatically; eliminating manual editing, reducing data error, and improving test quality to support critical patient needs.

### **Methods**

Since changing HLA Fusion to address the required format adjustment in an automated way will be challenging, we have identified a process to address it via the HistoTrac LABScreen custom import utility (OneLambda). A custom setting can be utilized to specifically address this need and promptly automate the DPB1 antibody data formatting adjustment during data import. Subsequently, the adjusted data is imported into HistoTrac, which can then be automatically uploaded to UNET via the UNOS API.

### Results

A custom setting in the LABScreen import utility has been employed (Figure 1) that successfully transforms the data to the UNOS acceptable format for all loci, including DPB1, regardless of how the call in the HLA Fusion software was made, thereby automating the data transfer between systems without the need for manual changes.

Figure 1. An optional custom setting in the HistoTrac LABScreen import utility to automatically adjust to the required DPB1 data format.

Import Settings	
Import General Information	Import CPRA
Import Test Results	Import PRA Result Code
Import Test Configuration Information	Import Specificity in Check Box for
Import Test Summary Information	🔲 Send Tests To Sign Off Queue a
Import Epitope Analysis	Import to Tests in the Sign Off Qu
Import Tail Analysis	Select Individual Samples to Imp
Import Manual Tail Analysis	Import DP Serology as 2-Field Mo
✓ Import Test Details	
<ul> <li>Import Test Details</li> <li>Split Bead Detail for Shared Beads</li> <li>Replace Notes on Import</li> </ul>	
<ul> <li>Import Test Details</li> <li>Split Bead Detail for Shared Beads</li> <li>Replace Notes on Import</li> <li>Import Tested Date From Fusion</li> <li>Use Build/Run Date vs Analysis</li> </ul>	Date for Tested Date
<ul> <li>Import Test Details</li> <li>Split Bead Detail for Shared Beads</li> <li>Replace Notes on Import</li> <li>Import Tested Date From Fusion</li> <li>Use Build/Run Date vs Analysis</li> <li>Valid Single Antigen Test:</li> </ul>	Date for Tested Date

Here we demonstrate how the automated adjustment works using a common example for DPB1. As shown in Figure 2, we assign DPB1 serological equivalents in HLA Fusion upon analysis.

When the assigned results are imported in HistoTrac via the LABScreen custom utility, with the new "Import DP Serology as 2-Field Molecular" option off, the specificities are displayed as equivalents (Figure 3).



Learn more at thermofisher.com/FURL

—		×	
t mport	<u> </u>	<u>D</u> K ose	
lar			



Figure 3. DPB1 serological equivalents in HistoTrac when the new "Import DP Serology as 2-Field Molecular" option is turned off.



When we turn on the new "Import DP Serology as 2-Field Molecular" option however, the specificities are displayed automatically as converted 2-Field molecular results (Figure 4).

Figure 4. DPB1 serological assignments automatically converted to molecular results in HistoTrac when the new "Import DP Serology as 2-Field Molecular" option is turned on. NOTE: The Specificity format displayed in the test code is listed as a per Locus string series representative of acceptable alleles/equivalents.



## **Conclusions**

Using custom software solutions to automate and improve processes has long been a vision of HistoTrac. The employed solution allows laboratories to meet the current nomenclature requirement of ASHI and UNOS, automating specificity calls without manual editing.

Additionally, it aids in minimizing the UNOS API and Waitlist data errors, ensures patient UNET profiles are accurate and up to date, saving precious time.

### **Trademarks/licensing**

© 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Science at a scan Scan the QR code on the right with your mobile device to download this and many more scientific posters.

# **Thermo Fisher** SCIENTIFIC

ugh the							-	0	×
Decult	operificity		Moderate Rick	I ow Pisk Antibody	Reportable Co	Notes	Save Data	Result	Comple
Result	DRw:3*01 DP:01:01 02:01 02 09:01 10:01 11:02	2:02 03:01 04:01 04:02 05:01 06:01 1 13:01 18:01 19:01 20:01 28:01		LOW RISK ARIBOUTY	Reportable Co	nutes		Kesun	
DR	wj DQB	DQA DPB DPA	View Cutoff ThresholdsApply Thresholds	Apply To Test	Spir Bea				
110 110	8						<u> </u>	_	
85 111 123 122	8 8 8	Specificity							
128 80 96	8 8 8	DRw:3*01	01 02:02 02:01 (	14-01-0	4.02	05-01-04			
89 110 107 104	6 6 6	09:01 10:01	11:01 13:01 18:0	)1 19:0	1 20:	01 28:01	L 1		
125 129 133	6 4 4								
129 118 111	4 4 4								
98 128 107 127	4 4 4								
107 113 114	4 4 4								
132 134 115	4								~
			Sample Comments		View Image View Details	Save And Close	<u>C</u> lose Wit	hout Sa et Test	wing
							<u>-1</u> 0		

