Comparison of two Elisa-based Angiotensin II -type-1 receptor (AT₁R) autoantibody detection assays using an antibody cross-reactivity check approach

Adak N. Karamafrooz¹, Robert Carroll², Dave Lowe¹, MichaelTrinh¹, Rui Pei¹, Julie Nguyen¹ 1. Research and Development, Thermo Fisher Scientific, West Hills, CA, United States. 2. Health Sciences, University of South Australia, Adelaide, SA, Australia.

Abstract

Elevation of autoantibodies against Angiotensin II Type 1 receptor (AT₁R) in human sera has been associated with several disorders, including pre-eclampsia and renalallograph rejection, therefore, an accurate and reliable measurement of those antibodies is crucial. The enzyme-linked immunosorbent assay (ELISA) is the primary method for measuring AT1R-specific antibodies (AT1R-Ab) levels. We have performed a cross-reactivity study comparing a commercial AT₁R antibody detection kit (EIA AT1RX, CellTrend) and a research-developed ELISA assay to demonstrate the accuracy and sensitivity of the ELISA as a detection platform for AT1R-specific antibodies.

Introduction

Rapid and accurate methods for the detection of human leukocyte antigens (HLA) and non-HLA antigens associated with the pathogenesis of hyperacute, acute, and chronic allograft rejection are clinically important for the successful management of long-term allograft transplantations. Among non-HLA antigens, antibodies against AT₁R have been extensively investigated and shown to be associated with vascular inflammation, leading to refractory vascular rejection and eventual allograft loss, enhancing the mortality risk in kidney transplant patients. In the current study, we verified the sensitivity of a commercial AT₁R-Ab detection kit (CellTrend) as well as a purified AT₁R mutant using AT₁R-specific antibodies that were eluted from live cells expressing AT₁R antigen. Furthermore, we discovered that certain reagents are incompatible with the commercial AT1R kit, resulting in a loss of detection signal.

Materials and methods

Transfection and expression of AT₁R: T-rex 293 adherent cells were maintained in DMEM high glucose supplemented with 10% FBS, 1% penicillin-streptomycin, and glutamax at 37°C and 5% CO2 in a humidified incubator. T-Rex 293 cells were stably transfected with AT₁R wild type (AT₁R-WT) or AT₁R constitutively active mutant using a Neon electroporation system. AT₁R surface expression was confirmed using flow cytometry and AT₁R-specific antibodies.

Purification of AT₁R using affinity chromatography: The purification was performed using affinity chromatography based on Wingler et al., 2019 with minor adjustments. AT₁R was eluted using 0.1 M glycine (pH 2.0-2.8) and neutralized with 1M Tris (pH 9.5). The eluted fraction was buffer exchanged with HEPES buffer (500 mM NaCl, 0.01% MNG, 0.004% (w/v) CHS, protease inhibitor, pH 7.5) and analyzed using BCA assay, western blot, and LC-MS mass spectrometry protein sequencing.

Preparation of ELISA tray for detection of anti-AT₁R autoantibodies: a) Stably transfected cells with AT₁R construct were seeded on Poly-D-Lysine-coated 96 well plates, induced 24 hours later using doxycycline, and fixed with %4 formaldehyde. b) The affinity-purified AT₁R antigen was supplemented with protease inhibitor and coated overnight at 4°C with shaking on the Maleimide-activated ELISA plates and then blocked with ELISA blocking buffer for 1 hour and washed three times before testing.

Antibody adsorption with crossmatch cells and elution (AXE procedure): T-rex293 cells expressing AT₁R were harvested 48 hours post-induction. Cell pellets were incubated with human serum positive for AT₁R-Ab for 30 minutes at 37°C with agitation. The supernatant was then retrieved, and the pellets were washed several times before antibody elution using elution buffer (Glycine buffer pH 3.0, MgCl₂) and neutralized using Tris buffer (pH 8.0). The serum-derived antibodies were tested on the AT₁R ELISA plates.

Adsorb out[™] bead (AOB) treatment: Adsorb Out[™] (One Lambda) microparticles are used to reduce the background reactivity due to nonspecific binding. Briefly, test serum and AOB were dispensed in varying dilutions, and incubated for 30 minutes at room temperature with agitation. Following centrifugation, the supernatant was transferred into a new tube without disturbing the pelleted beads. The sera were then transferred to the ELISA tray for further processing.

Results





Figure 2: effect of BSA on the CellTrend AT₁R kit's detection signal. a). BSA, as the major component of the AOB buffer, impairs the detection signal in the positive control of the AT₁R Kit, while other components do not interfere with the signal detection. **b**): The presence of BSA, even at very low concentrations, affects the signal detection on the CT Elisa plate. Samples are titrated with varying BSA concentrations prior to testing on the ELISA tray.



Figure 3: Cell-based ELISA tray: Standards and controls from CellTrend AT₁R kit as well as patients' sera reacted to the ELISA tray Coated with T-Rex 293 Cells expressing AT₁R active mutant, confirming the specificity of AT₁R antigen present on the cell surface.





Figure 1: Effect of serially diluted Adsorb Out bead (AOB) on standards & controls of the CellTrend AT₁R-Ab Assay kit (lot 64). Serially diluted AOB was added to the kit's standard sample, positive and negative controls, or diluent buffer (DB). When used at high concentrations, AOB would interfere with the detection signal. At higher dilutions (10^{-7}) , the inhibitory effect is diminished, and the signal is restored.

Figure 4: Purified AT1R-coated ELISA tray. Affinity-purified AT₁R consecutively active mutant is coated on an ELISA tray. *CellTrend* standards (4-a) as well as human serum positive and negative for AT₁R-Ab (4-b) are tested on the AT1R-coated ELISA tray.



CellTrend standard (U/ml) and controls

Figure 5: Cross-reaction of the eluted AT₁R-specific antibodies with the ELISA trays. Using AXE method (Ref.1), eluted antibodies from patient serum positive for AT₁R, which had been captured on AT₁R-expressing cells, cross-reacted with both the AT₁R ELISA tray from the CellTrend kit (5-a) and the purified AT₁R coated Elisa plate (5-b).





Thermo Fisher S C I E N T I F I C

Conclusions

- There have been studies that show pretreatment of the samples with Adsorb out[™] bead prior to performing the ELISA assay can lead to a decrease or complete loss of detection signals (reviewed in Ref. 2). Here we demonstrate that bovine serum albumin (BSA) present in the AOB buffer is the primary element that impedes the detection signal in the commercial AT₁R ELISA kit.
- Sensitivity and specificity of ELISA as a reliable platform are confirmed through the Adsorption with Crossmatch cells and Elution (AXE) experiment. This experiment defines only the antibodies bound to the specific antigens on the cell surface. By transferring the eluted antibodies to the AT₁R-Ab detection ELISA tray, we successfully detected those antibodies.

References:

- Leahy, R., Sullivan, L., & Carroll, R. (2023). P403 The adsorption crossmatch cells and elution (AXE) technique to identify true HLA-specific antibodies. Human Immunology, 84, 122.
- 2. Xu, Q., Johnson, K. P., Hardiman, M., Helmick, D., & Zeevi, A. (2023). Loss of anti-AT1R reactivity in ELISA post-adsorption - False reactivity or interference in the assay? Human Immunology, 84(4), 286–289.
- 3. ThermoFisher Sci, (n.d.). Deciphering the Loss of Anti-Angiotensin Type 1 Receptor Reactivity in Enzyme-Linked Immunosorbent Assay Post-Adsorption. P519 Human Immunology, Volume 84, Supplement, September 2023, p142-143.

Acknowledgements:

We thank Professor Robert Carroll from the University of South Australia for supplying the human sera.

Trademarks/licensing

© 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.