An HLA genotype-informed dPCR assay for monitoring donor-derived cfDNA

¹Bram Luiken, ¹Lara van der Klugt, ¹Bo Rood, ²Monica Sorbini, ²Silvia Deaglio, ²Tiziana Vaisitti, ¹Bart Valkenburg, ¹Freek Manders, ¹Eleni Draina, ¹Sake van Wageningen

¹GenDx, Utrecht, The Netherlands

²Department of Medical Sciences, University of Turin, Italy

Selection assay Pre-amp multiplex Digital PCR Data analys



-- Homozygous Donor Mismatched Allele -- Heterozygous Donor Mismatched Allele

Results:

The current protocol shows a correct trending of the FAM signal (R2>0.99 for heterozygous samples from 0.2% to 100% and R2>0.98 for homozygous from 0.2% up to 50% (Figure 3a&b, Figure 4), but requires further loading optimization for the homozygous FAM and for the HEX channel in general.

Conclusion:

Sensitive and specific results can be achieved using SOtrack. The use of HLA markers allows for a rapid test with a turnaround under 3.5 hrs without the need for a pre-typing step.

Aim:

Donor-derived cell-free DNA (cfDNA) is being evaluated as a biomarker for monitoring graft injury in solid organ transplantation. Monitoring cfDNA levels in peripheral blood could be used to confirm or substitute biopsies. The assay presented here, SOtrackTM, uses markers based on HLA typing in combination with dPCR. we test the performance of HLA-DRB1 and HLA-DQB1 allele group-specific SOtrack markers using the Quantstudio Absolute Q dPCR system on artificial chimeric cfDNA samples.

