

# The antimicrobial activity of silver-containing dressings against surface-associated microorganisms in a dual-species community

Kate Meredith, PhD<sup>1</sup>, Matilda Coleborn, MSc<sup>1</sup>, Lucy Forbes, BSc<sup>1</sup>, Daniel Metcalf, PhD<sup>1</sup>

<sup>1</sup>Convatec Ltd, Deeside, UK

## Introduction

- Microorganisms found in hard-to-heal wounds are invariably present in a surface-attached phenotype<sup>1</sup>
- These surface- and self-attached microbial communities usually contain multiple species<sup>2</sup>
- This polymicrobial nature of wound bioburden increases virulence and tolerance to antimicrobial agents<sup>3,4</sup>
- We evaluated antimicrobial dressings containing silver against surface-associated microorganisms in a dual-species community grown in a stringent simulated wound model

## STUDY OBJECTIVE

The aim of this *in vitro* study was to evaluate the effectiveness of different antimicrobial dressings against dual species communities grown in a simulated wound model

## Methods

### Test dressings

- Carboxymethylcellulose dressing containing ionic silver, ethylenediaminetetraacetic acid (EDTA), and benzethonium chloride (BEC) ('**CISEB**')
- Non-adherent polyethylene mesh with polyester core dressing containing silver oxysalts ('**PPSO**')
- Cellulose ethyl sulphonate fibre dressing containing ionic silver ('**CESIS**')
- Polyacrylate (polyabsorbent) fiber dressing with acrylic core containing silver sulphate ('**PSS**')
- Secondary film cover dressing

### Challenge organisms


- Klebsiella pneumoniae* (NCTC 9156)
- Community-associated Methicillin-resistant *Staphylococcus aureus* (CA-MRSA; ATCC® BAA™-1556, clone of USA300)

### Protocol


- Challenge organisms were grown on gauze for 24 hours in a nutritious broth and fetal bovine serum (FBS) mixture with an initial inoculum number of  $\sim 1 \times 10^6$  colony forming units (CFU)/mL in total
- Following incubation, the gauze-attached dual-species community was rinsed twice to remove planktonic bacteria and placed on the central contact plate within a simulated wound assembly (SWA), and total viable counts were performed to establish initial bacteria numbers
- The SWA consisted of a bovine leather-covered Perspex plate (simulating peri-wound skin), surrounding a central insert of a 55 mm diameter agar contact plate (simulating a moist wound bed with a reservoir of isotonic nutrients), supporting the dual-species communities
- The simulated wound area was covered with the test primary dressing (10×10 cm samples) and hydrated with 8 mL simulated wound fluid (50/50 Maximum Recovery Diluent: FBS)
- The SWA was incubated at  $35 \pm 3$  °C for 6, 24, 48, and 72 hours
- At each timepoint, enumeration of surviving bacteria was performed, with n=3 for each test dressing and n=1 for the equivalent control; whereby no dressing was applied to the surface attached dual species community

## Results

- CISEB** dressing was observed to reduce the dual-species community within 24 hours with a  $>4 \log_{10}$  kill observed for *K. pneumoniae* and a  $>6 \log_{10}$  for MRSA from an initial challenge of  $\sim 1 \times 10^8$  CFU/mL for each separate challenge organisms present
- This kill rate was sustained for the duration of the challenge period, with **CISEB** dressing reducing the dual-species community to non-detectable levels ( $<30$  CFU per test) by 72 hours for *K. pneumoniae* and by 48 hours for MRSA
- PPSO** dressing demonstrated an initial reduction at 6 hours of  $\sim 2 \log_{10}$  for both *K. pneumoniae* and MRSA. However, this dressing was unable to further reduce the population after this time, with bacterial numbers at subsequent timepoints being similar to or greater than those at 6 hours for both bacterial species
- CESIS** dressing exhibited a slight, gradual reduction in bacterial numbers for the dual-species community over the course of the test period, reducing each challenge organism by  $\sim 2.5 \log_{10}$  in 72 hours
- PSS** dressing was shown to have little to no impact on the dual-species community, with levels remaining similar or greater than those recovered before dressing application, the initial growth (the  $T_{0hr}$  count)
- The no-dressing **control** demonstrated that the dual-species community remained viable throughout the test period and that species population proportionality was maintained



**Figure 1.** Antimicrobial activity of test dressings against a surface-attached dual-species community (total population) of *K. pneumoniae* and CA-MRSA over a 72-hour test period.



**Figure 2.** Antimicrobial activity of test dressings against a surface-attached dual-species community (challenge organisms shown individually) of *K. pneumoniae* and CA-MRSA over a 72-hour test period. The data has been normalized by subtracting the limit of detection (30 CFU).

## Discussion

- The surface-attached dual-species community considers the polymicrobial nature of microbial communities found clinically, and this study showed that not all silver-containing dressings are equally effective against these complex communities
- It is important that wound dressings contain proven additional components, that can aid in the breakdown of these surface-attached communities, as well as an antiseptic agent, to ensure the effective reduction of polymicrobial surface-attached communities.

## CONCLUSION

Using a surface-attached dual-species community in a stringent *in vitro* test, it was demonstrated that not all silver-containing dressings can reduce these types of microbial communities. Only CISEB showed enhanced activity, which may be attributed to the additional agents contained within the dressing – EDTA and BEC – in addition to ionic silver

## References & Footnotes

- Malone M et al. *J Wound Care* 2017;26(1):20–25.
  - Wolcott R et al. *Clin Microbiol Infect* 2013;19(2):107–112.
  - Anju VT et al. *Antibiotics (Basel)* 2022;11(12):1731.
  - Orazi G & O'Toole GA. *J Bacteriol* 2019;202(1):e00530-19.
  - Bowler P & Parsons D. *Wound Medicine* 2016;14:6–11.
- CISEB:** Aquacel® Ag+ Extra™ (Convatec, UK)  
**PPSO:** KerraContact® Ag (Crawford Healthcare, UK)  
**CESIS:** Durafiber Ag (Smith & Nephew, UK)  
**PSS:** UrgoClean Ag (Urigo Medical, UK)  
Secondary film cover dressing: Tegaderm™ (3M Health Care, UK)