

DETECTION OF PSEUDOMONAS AERUGINOSA CONCENTRATION VIA SYNTHETIC BIOSENSOR FOR QUORUM SENSING AUTOINDUCERS




Skylar Leslie¹, Jordon Gilmore¹, David Karig¹
(1) Department of Bioengineering, Clemson University



INTRODUCTION


Pseudomonas aeruginosa is an opportunistic pathogen commonly found across various environments and hosts (1, 2). *P. aeruginosa* prevalence (32,000 infections in hospitalized patients, 2,700 deaths in 2017 (3)) presents a severe threat with multiple drug resistance (MDR) (1). Contributing to MDR, *P. aeruginosa* coordinates biofilm formation through quorum sensing (QS), in which diffusible molecules are secreted and sensed to detect local population density (2,4). This work focuses on detecting *P. aeruginosa* QS signaling molecules to inform infection diagnosis and treatment. Specifically, we detect N-(3-Oxydodecanoyl)-L-homoserine lactone (3OC12HSL) from the *las* QS pathway, which coordinates virulence. A portion of the *P. aeruginosa* pathway can be seen in the figure below (5).



3OC12HSL binds the protein LasR, forming a complex that activates Las QS-controlled promoters. We employ a synthetic plasmid that codes for the expression of LasR from a constitutive promoter and the expression of a colorimetric (blue) reporter, *amiCP*, from a QS-controlled promoter. The biosensor host is nonpathogenic *Escherichia coli*, offering a sensitive, specific, and low-cost detection method with visual readouts interpretable without additional equipment.

METHODS

Transformation



- 3uL of GA10-10 *amiCP* vector (seen above) was added to the DH5a *E. coli* cells and incubated on ice for 20 minutes.
- Cells were then heat shocked at 42°C for 30 seconds. The cells were then placed back on the ice for 2 minutes.
- 400uL of SOC broth was added, and the solution was transferred to a 14mL tube and then placed in a 37°C incubator for 1 hour.
- Cells were then plated on LB + kanamycin agar and grown overnight.
- For later testing, competent cells were frozen in a 30% glycerol stock solution.

Testing

Colorimetric Testing:

- Transformed DH5a+*amiCP E. coli* cells were then cultured in LB+kanamycin for 18 hours.
- Conical tubes (n=12) were prepared as follows: 4mL with 10uM exogenous 3OC12HSL added to tube 1; 3mL with 1 mL serially diluted aliquots added to tubes 2-12.
- Tubes were plated in 48-well plates, placed in a 37°C plate reader, and read at 588nm, 600nm, and 700nm for 24 hours.
- In the presence of 3OC12HSL, DH5a+*amiCP E. coli* cells turn blue, providing easy visual readout without equipment.

Media Testing:

- A *P. aeruginosa* colony was inoculated in 3mL of PBS and vortexed.
- Seven different culture types were prepared: 25-5% LB + 75-95% artificial wound fluid exudate (AWFE), 3mL each.
- 300uL of *P. aeruginosa* in PBS was mixed into each solution.
- 200uL of each solution was added to a 96-well plate in four replicates. The plate was read every 30 minutes at 600nm, 37°C 200RPM for 24 hours.
- After 24 hours, a Crystal Violet assay was performed to measure biofilm formation at 590nm.

Sol-gel Preparation

- 5mL of 0.4M Sodium Silicate (SiNa) and 5mL of 6.5M of LUDOX were added to a 250mL beaker and mixed at 200 RPM for 30 seconds.
- 5mL of overnight DH5a+*amiCP E. coli* and 10mL of Phosphate buffer solution (PBS) were mixed for 1 minute in the beaker.
- Initial pH was measured (around 10.8), and 1M HCl was added to the solution until pH was between 7.0 and 7.3.
- 1mL aliquots were placed in a 24-well plate and allowed to gel.
- Gelation of solution occurs within 12 minutes of the addition of HCL.

RESULTS




Figure 1: Detection of 3OC12HSL. Blocks indicate different concentrations of exogenous 3OC12HSL with dilutions starting at A1:A3 and ending at F4:F6. Bacterial control is seen at A7:C7 and media control at D8:F8.




Figure 2a: Biosensor growth curves when exposed to varying concentrations of exogenous 3OC12HSL.




Figure 2b: Response curve showing the concentrations of exogenous 3OC12HSL normalized to bacterial control over time.




Figure 4a: *P. aeruginosa* growth curves based on varying LB concentrations.




Figure 4b: Biofilm analysis of different media concentrations. $P=0.05$ for all statistical tests.



Figure 3: Encapsulation of *amiCP* biosensor. The reaction of SiNa and LUDOX with HCl allows for room-temperature encapsulation of bacteria.

$$\begin{aligned} \# \text{ of molecules} & \times \frac{1}{6.022 \times 10^{23}} = 1.66 \times 10^{-24} \text{ moles} \\ (X) \times 1.66 \times 10^{-24} \text{ moles} & \times \frac{1}{6 \times 10^{16} \text{ L}} = (X) \times 2.77 \times 10^{-9} \text{ M} \\ 0.000152588 \mu\text{M} & \times \frac{1 \times 10^{-6} \text{ M}}{1 \mu\text{M}} = 1.52588 \times 10^{-10} \text{ M} \\ (X) \times 2.77 \times 10^{-9} \text{ M} & = 1.52588 \times 10^{-10} \text{ M} \\ X & = 0.055 \text{ molecules/cell} \end{aligned}$$

Equation 1: Equation used to determine the number of 3OC12HSL molecules per biosensor cell at 0.15nM concentration.

DISCUSSION AND CONCLUSION

- DH5a+*amiCP E. coli* cells could detect exogenous 3OC12HSL in concentrations as low as 0.15nM.
- Growth curves and response curves can be seen in **Figures 2a** and **2b**.
 - 3OC12HSL may have some effect on biosensor growth.
 - We can see the correlation between **Figures 1 & 2**.
- Calculations in **Equation 1** indicate that our biosensor can detect as low as 0.055 molecules of 3OC12HSL/cell.
- DH5a *amiCP E. coli* could be encapsulated in sol-gel, as seen in **Figure 3**.
- *P. aeruginosa* growth and biofilm development depend on the LB concentration in liquid media (**Fig. 4a** & **Fig. 4b**).

FUTURE DIRECTION

- Usage of mass spectrometry may allow us to quantify the amount of 3OC12HSL molecules present in relation to bacterial density.
- Cell viability testing and biosensor validation within the sol-gel solution will allow us to determine if sol-gel is a valid method of detecting 3OC12HSL molecules.
- Testing the biosensor against polymicrobial communities will further validate its ability to detect 3OC12HSL molecules.

REFERENCES

1. 2019. *Pseudomonas aeruginosa* Infection | HAI | CDC.
2. Thi MTT, Wibowo D, Rehm BHA. 2020. *Pseudomonas aeruginosa* Biofilms. Int J Mol Sci 21:8671.
3. CDC. 2022. The biggest antibiotic-resistant threats in the U.S. Centers for Disease Control and Prevention.
4. G. Glansdorp F. L. Thomas G, K. Lee J, M. Dutton J, C. Salmond GP, Welch M, R. Spring D. 2004. Synthesis and stability of small molecule probes for *Pseudomonas aeruginosa* quorum sensing modulation. Organic & Biomolecular Chemistry 2:3329–3336.
5. Soukariéh F, Williams P, Stocks MJ, Cámará M. 2018. *Pseudomonas aeruginosa* Quorum Sensing Systems as Drug Discovery Targets: Current Position and Future Perspectives. J Med Chem 61:10385–10402.

Acknowledgments: The authors want to acknowledge Craig Miller, Diptee Chaulagain, and the SCINBRE DRP Award for their assistance with previous work and funding opportunities.